



## Defibrillator Safety & Performance Testing

DUT (Device Under Test): \_\_\_\_\_

Objective: To collaboratively test and verify the safety and performance of a Defibrillator in various modes (Manual, Synchronized Cardioversion, Pacing) using a defibrillator analyzer, adhering to electrical safety and procedural standards.

Pre-settings: Ensure the analyzer and DUT are prepared and set ready for testing. Confirm electrical safety, proper lead connections, and correct analyzer mode selection.

### Part A: Electrical Safety Tests (IEC 60601-1)

DUT Details: Model: \_\_\_\_\_ Serial No.: \_\_\_\_\_

Test Leads: Patient Cable Type: \_\_\_\_\_ Test Load: 50Ω

Test	Parameter / Mode on DUT	Measurement from Impulse D6000	Expected Range / Limit	Pass/Fail
Energy Accuracy Test	Manual Mode: Set: 50J, 100J, 200J Charge & Deliver	Measured Energy (J): 50J: _____ J 100J: _____ J 200J: _____ J	Tolerance: ±15% or ±4J (whichever is greater)	[ ] Pass [ ] Fail
Synchronized Cardioversion	Sync Mode: Set: 50J Sync to R-wave & Deliver	Energy Delivered: _____ J Sync Delay: _____ ms	Energy: ±15% Delay: < 60 ms from R-wave	[ ] Pass [ ] Fail

Team Observations & Troubleshooting Notes:

---

---

---

---

---

---



Learn Beyond

(Autonomous, Affiliated to Anna University)

**KPR Institute of  
Engineering and  
Technology**

**KPR Center for Biomedical  
Equipment Testing and Calibration**

## Charge Time Test

*Time taken to charge to selected energy*

Energy Setting	Maximum Allowable Charge Time	Measured Charge Time	Pass/Fail
50 J	15 seconds	_____ sec	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
100 J	15 seconds	_____ sec	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
150 J	15 seconds	_____ sec	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
200 J	20 seconds	_____ sec	<input type="checkbox"/> Pass <input type="checkbox"/> Fail

## Custom ECG Waveform Generation for performance test

Test	Custom Waveform Parameters (Generated in Impulse D6000)	Expected Result	Pass/Fail
1. VF Detection - Fine Fibrillation	Amplitude: <b>0.2-0.5 mV</b>	Shock Advised = <b>Yes</b> Time < 15 sec	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
2. VF Detection - Coarse Fibrillation	Amplitude: <b>1.0-2.0 mV</b>	Shock Advised = <b>Yes</b> Appropriate energy	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
3. VT with Pulse Simulation	Rate: <b>180 bpm</b>	Manufacturer dependent (Often "No Shock")	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
4. Asystole Detection	Flatline: <b>&lt; 0.1 mV</b>	"No Shock Advised" CPR Prompt = Yes	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
5. Sinus Rhythm + Artifact	NSR 80 bpm + <b>Motion Artifact</b>	Should NOT advise shock	<input type="checkbox"/> Pass <input type="checkbox"/> Fail



Learn Beyond

(Autonomous, Affiliated to Anna University)

**KPR Institute of  
Engineering and  
Technology**

**KPR Center for Biomedical  
Equipment Testing and Calibration**

## AED Algorithm Performance Scoring

\*Rate the AED's performance on each custom waveform (1=Poor, 5=Excellent)\*

Waveform Type	Detection Accuracy	Analysis Speed	Overall Score
Fine VF	[1] [2] [3] [4] [5]	[1] [2] [3] [4] [5]	_____ / 10
Coarse VF	[1] [2] [3] [4] [5]	[1] [2] [3] [4] [5]	_____ / 10
Any other arrhythmia	[1] [2] [3] [4] [5]	[1] [2] [3] [4] [5]	_____ / 10
Noisy Signal	[1] [2] [3] [4] [5]	[1] [2] [3] [4] [5]	_____ / 10
<b>TOTAL</b>			_____ / 40

## Post-Test Analysis & Reporting

Data Analysis:

Did any test point fail? If yes, what could be the potential causes?

---

---

---

---

Error Sources:

List three potential sources of error in defibrillator testing:

---

---

---

Clinical Significance:

Why is energy accuracy ( $\pm 15\%$ ) critical in defibrillation?

---

---

---



Learn Beyond

(Autonomous, Affiliated to Anna University)

**KPR Institute of  
Engineering and  
Technology**

**KPR Center for Biomedical  
Equipment Testing and Calibration**

### Team Handover & Summary:

- Defibrillator model tested: \_\_\_\_\_
- Overall result (All Pass / Fail at X test): \_\_\_\_\_
- One key challenge faced and how it was resolved: \_\_\_\_\_